




GL 60 mg/dL
& mmol/L

beurer
medical

GL 60 mg/dL & mmol/L Blood glucose monitor



Specification of a target range for visual classification of the measurement result



Slim, modern design with integrated USB connector



Compatible with the free "beurer HealthManager Pro" app



Starter set: Measuring device GL 60, lancing device LD 04, 10 lancet needles, 10 test strips, case, instructions for use, short instruction, blood glucose diary

Backlit display in landscape format

Dextrose as well as ketone warning on display

Measured value memory for 900 measurements with date, time and measurement marking

Easy measurement marking for optimized diabetes management: before eating, after eating, general marking

Calculating the average value for 7, 14, 30 or 90 days for
- all measurements
- marked measurements

Wide hematocrit range from 0-70%

Blood sample: fresh whole blood (capillary) or venous whole blood

Required volume of blood: 0.5 microlitres

Blood volume check: measurement only possible with sufficient sample of blood

Duration of measurement: approx. 5 seconds

Measurement range: 10 - 630 mg/dL / 0.56 - 35.0 mmol/L

Automatic switch-off: 2 minutes after last operation

Suitable for self-testing

Code-free

Plasma-calibrated

Integrated USB plug

Data transfer via *Bluetooth*[®] and USB

2 x 3 V CR 2032 Button cell batteries included

Product dimensions: 112.6 mm x 29.9 mm x 14.7 mm

Product weight: approx. 40 g

Fully fulfils EN ISO 15197:2015





GL 60 mg/dL
& mmol/L

beurer
medical

Storage/transport:	Temperature:	Rel. humidity:
Device:	-20 °C bis +50 °C	10 – 95 %
Test strips:	+2 °C bis +32 °C	<90 %
Operation:	+10 °C bis +40 °C	30 – 85 %

Medical device, 5-year guarantee

Sales unit: 5 / Shipping carton: 4 x 5

Item no.: 463.55 (mg/dL)	463.57 (mmol/L)
EAN no.: 4211125463557 (mg/dL)	4211125463571 (mmol/L)
PZN no.: 16382624 (mg/dL)	16382630 (mmol/L)
HMV no.: XX.XX.XX.XXXX (mg/dL)	XX.XX.XX.XXXX (mmol/L)



List of system requirements and compatible devices

GL60_2023-06-12_01_DS_EN_BEU Subject to errors and changes

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Beurer GmbH is under license. Other trademarks and trade names are those of their respective owners.

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries.

Google Play and the Google Play logo are trademarks of Google LLC. Android is a trademark of Google LLC.

*This product satisfies the requirements of the applicable European directives.